

Nasal Suction Therapy at Home

Training and Safety Guidelines

For your safety

Following these guidelines will help reduce potential risks. See the manufacturer website for more information about your equipment.

- Plug the three-prong power cord only into outlets that are properly grounded. An electrician can show you the grounded outlets in your home.
- Place the machine on a hard, flat surface.
- Clean all equipment regularly.
- Check filters regularly.
- Wash hands regularly.

How to use the suction machine

Setup

- 1 Place the suction machine on a stable surface and plug it in.
- 2 Ensure the canister is attached properly to the unit as shown above.
- 3 Connect the tubing to the outlet port on the lid of the collection bottle.
- 4 Attach the suction device to the other end of the tubing.



Therapy

- 1 Wash your hands.
- 2 Turn on the machine.
- 3 Turn the pressure adjustment knob to the setting directed by your healthcare provider.
- 4 Gently insert the catheter into each nostril and suction for no more than 10 seconds.
- 5 After each nostril, suction water through the catheter until the catheter and tubing are clear.
- 6 Turn off the machine.
- 7 Wash your hands.

Clean Regularly

Follow these guidelines to keep your suction equipment clean and help prevent infection.

Daily

Empty the collection bottle at least once each day — flush fluids down the toilet. Wash the bottle, lid, tubing, and catheter in warm soapy water and rinse well. Allow to air dry.

Weekly

- Wipe the outside of the machine with a damp cloth.
- After washing the suction bottle, lid, tubing, and catheter with soap and water, soak them in a water and vinegar solution (3 parts water to 1 part vinegar) for at least 20 minutes. Rinse well and air dry.

Every Few Weeks

- Regularly check the filters and replace them as needed, at least monthly.
- Replace tubing at least monthly or when it can no longer be cleared adequately.
- Replace catheters at least monthly.

Follow common troubleshooting tips

If your suction machine is not working, check these things:

- Is the machine plugged in?
- Is the power switch turned on?
- Does the wall outlet the machine is plugged into have power?
- Is the collection bottle full?
- Is the tubing securely fastened?
- Is the collection bottle top securely fastened?

Call Intermountain Homecare & Hospice if...

You experience any problems with your equipment.

Notes

Content created and approved by Intermountain Homecare & Hospice

Intermountain Healthcare complies with applicable federal civil rights laws and does not discriminate on the basis of race, color, national origin, age, disability, or sex. Se proveen servicios de interpretación gratis. Hable con un empleado para solicitarlo. 我們將根據您的需求提供免費的口譯服務。請找尋工作人員協助。