



# USING A BROWN BAG WITH HANDLES FOR N95 RE-USE AND DISINFECTION

05/27/20 v2

1


Unfold and flatten the side creases of a brown bag with handles so it stays open.



Opening the bag better supports hanging the N95 by handles during re-use throughout the day.

2

When you need to remove your N95, you can hang it inside the bag by looping the strap around the handle.



This will limit your potential contact with the face piece of your N95 as you reuse it throughout your shift.

3

At the end of your shift, place the labeled N95 at the bottom of the bag with handles.



Caregivers who use N95 masks but disinfection is unavailable to them should also use these brown bags with handles.

Fold the bag at the bottom of the handles.



Then place the labeled bag in the tan tote for reprocessing.

