## **High Risk COVID-19 Isolation:**

## **Contact/Airborne PAPR**



March 7, 2020

This type of isolation should be used for patients who are classified as 'High Risk' and require testing.

## **Donning:**

- 1. Donning takes place outside the patient care area or in an anteroom.
- 2. Before donning, perform hand hygiene.
- 3. PAPR
  - a. Inspect motor blower, battery level, breathing tube, filter, and headcover.
  - b. Turn on blower.
  - c. Perform flow check with flow indicator.
  - d. Attach breathing tube to blower and headcover.
  - e. Disinfect inside of headcover with wipe (purple, gray, or orange top wipe).
  - f. Fasten PAPR blower around the waist.
  - g. Place headcover over your head and ensure it is sealed underneath chin and face.
- 4. Gown
  - a. Don isolation gown.
  - b. Fasten behind neck and waist.
- 5. Gloves
  - a. Don gloves.
  - b. Ensure glove covers cuff of isolation gown.

## **Doffing:**

- 1. Doffing of gown and gloves takes place just inside the patient room at the door, before leaving the care area. PAPRs are removed outside of the patient care area.
- 2. Gown and Gloves
  - a. External gown and gloves may be contaminated. Avoid contamination of self, others, and environment when removing. If grossly contaminated, wipe PPE with disinfectant wipe.
  - a. Unfasten gown ties or break by pulling carefully at the front of the gown.
  - b. Carefully, peel off gown and gloves from the inside out, roll up and discard.
  - c. Perform Hand Hygiene.
  - d. Exit patient care area.
- 3. PAPR
  - a. Perform Hand Hygiene.
  - b. Don clean gloves.
  - c. Carefully remove headcover by pulling away from face at chin level and set on cleanable surface.
  - d. Disinfect using purple top, gray top wipes, or orange top wipes.
    - i. Ensure the wet/contact times are met according to the wipes used:
      - 1. purple top (2 min) gray top (3 min), orange top (4 min)
  - e. Wipe down inside of headcover.
  - f. With a new wipe, clean outside of headcover.
  - g. With a new wipe, clean breathing tube.
  - h. With a new wipe, clean belt and blower.
  - i. Place PAPR in storage area.
  - j. Clean work surface and wipe dispenser with disinfectant wipe.
  - k. Remove gloves.
  - I. Perform Hand hygiene.