


CADD Infusion Pump for Continuous Chemo

Your doctor has prescribed a medication that you will give yourself using a CADD infusion [in-fyoo-zhun] pump. You have received detailed instructions with your pump. Use this sheet as a quick reminder for some of the things you learned in your instructions.

About your CADD infusion pump


Using the keypad

Keypad buttons on your pump and how they are used:



Home screen information

What you will see on the screen:



- If the IV bag begins to leak, place it in a plastic bag and call.
- Pump should not get wet. Place backpack outside of shower on hook during shower/bath time.
- Do not store sharp items in pump backpack.
- For "High Pressure" alarms: Check for kinks in tubing.
- Keep fresh AA batteries on hand with pump.

If you need to stop the pump:

You may need to stop the pump to change the batteries, if you're having trouble fixing the alarms with the help screen, if you think you're having a reaction to the medication, or there is a problem with the line. Follow these steps to stop the pump:

- 1 Press the **stop/start** key.
- 2 Press the **on/off** button on the side to power off the pump.
- 3 Select Yes when the **Power down?** message appears. The pump will turn off.
- 4 Clean your hands.
- 5 Disconnect your IV if your nurse has instructed you to do this yourself.

Your pump's batteries

Your pump comes with a **rechargeable battery pack** that you will use most of the time.

Recharge the pump regularly. When the display reads **Battery Low, Battery Depleted, or Rechargeable battery reached end of use**, plug the pump into the AC adapter. You may use the pump while it is recharging.

If you are away from a power source and can't plug into the AC adapter, you may use normal AA batteries (but not rechargeable or heavy duty batteries). Change all 4 batteries at the same time. Do not mix old and new batteries.

To change batteries:

- 1 Turn off the pump.
- 2 Using your fingers or a coin, rotate the knob on the battery door on top of the pump.
- 3 Remove the battery pack and insert 4 new AA batteries. Replace the battery pack and battery door afterward.



ALARMS: Your pump has several alarms.

- **Acknowledge** clears the alarm.
- **Silence** shows the alarm screen but silences the alarm for 2 minutes. This continues until the alarm is acknowledged or resolved.
- If help screens are available for the alarm, select **Help** and follow the instructions.
- Call Infusion Services or your healthcare provider if alarms continue.

Need help or have questions?

- Nurse: _____
- Phone number: _____
- Between 8:00 a.m. and 5:00 p.m., contact Infusion Services. See the table below for your location's phone number.
- Between 5:00 p.m. and 8:00 a.m. (after hours), contact your physician's office on-call provider for instructions.

Alta View Infusion	801-501-2436
Cedar City Cancer Infusion	435-868-5690
Intermountain Medical Cancer Infusion Center	801-507-3980
Layton Infusion	801-543-6280
Logan Cancer Center	435-716-6400
Loveland (LDS) Infusion Clinic	801-408-3540
McKay-Dee Cancer Infusion Center	801-387-7162
Primary Children's Cancer Center	801-662-4700
St. George Cancer Center	435-251-4900
Utah Valley Infusion Center	801-357-7341